User's Instructions Gebrauchsanleitung Manuel de l'utilisateur
Gebruiksaanwijzing Manuale utente Bruksanvisningar


INNOLUX BRIGHT

Innojok Oy
Sirrikuja 3 L
00940 Helsinki
Finland
p. +358(0)9 4789 2200
f. +358(0)9 04789 2220
www.innojok.fi

Innojok Oy Sirrikuja 3 L 00940 Helsinki Finland p. +358(0)9 4789 2200 f. +358(0)9 04789 2220 www.innojok.fi

INNOLUX BRIGHT

English.....	2
Deutsch.....	4
Français.....	6
Nederlands.....	8
Italiano	10
Svenska.....	12
Suomi.....	14
Eesti keeles.....	16
Русский.....	18
Latviešu	20
Lietuvių k.	22

English.....	2	Deutsh.....	4	Français.....	6
Nederlands.....	8	Italiano	10	Svenska.....	12
Suomi.....	14	Eesti keeles.....	16	Русский.....	18
.....	20	Lietuvių k.	22	Latviešu	

V. 13A

Congratulations for purchasing INNOLUX bright light therapy lamp. Before using the lamp, it is good to know what effects bright light has on your body and how to use bright light therapy effectively. Many people suffer from lethargy and depression during the dark winter season. This condition is known as SAD (Seasonal Affective Disorder). Scientific research has shown that as the amount of natural light decreases, the body may physiologically respond to these "unfavourable living conditions" by adjusting certain hormone levels and biological activities. SUB SAD or the Winter Blues is a milder form of SAD affecting a considerably larger number of people.

The effects of these adjustments may include:

- Intense fatigue, exhaustion and various sleep disorders
- Increased appetite leading to weight gain
- General passivity and depression
- Decreased work efficiency and loss of vitality

On a sunny day the strength of light is approximately 50 000 – 100 000 lux, whereas indoors it is as little as 100-200 lux. An effective bright light therapy session requires the eyes to be exposed to at least 2 500 lux. Via the retina and the optic nerves the light affects the brain's light sensitive region (suprachiasmatic cellkernel), which adjusts the body's physiological cycles in accordance with the rhythm of day and night. The light reduces the amount of the sleep hormone melatonin, secreted by the pineal gland. Melatonin is a key factor in the regulation of alertness and sleep cycles, but also in the prevention of SAD. Bright light therapy helps you alleviate these problems and regain vitality. If the symptoms are severe, they should be treated under medical supervision.

THErapy INSTRUCTIONS

Start bright light therapy before any symptoms of SAD appear, preferably soon after nights start drawing in.

Take a bright light session 3-5 times a week throughout the winter. If the symptoms are not relieved or they reappear, try taking daily sessions for 1-2 weeks. You may even begin with daily sessions or repeat this intensity whenever you feel you need it.

One session lasts for 0,5-2 hours. Longer than two hours does not intensify the effect but there's no harm of longer sessions. It is advisable to take the session between 6.00 a.m. and 6.00 p.m. but it is most effective in the morning between 6 and 8. A session in the evening may disturb sleep. If you tend to have evening sleepiness (common to elderly people), you might want to try taking the sessions earlier in the evening in order to stay up longer. You can also divide the session into several shorter periods during the same day if it suits your schedule better.

There is no need to look straight into the light source, but you should keep it within your visual field during the whole session. You can read, write or do other things in the meantime. The effectiveness of bright light therapy correlates with the distance between your face and the light source. Make sure that your face is not further from the lamp than the announced therapy distance on the model instructions. You can shorten the required session time by sitting closer to the lamp.

Bright light therapy also helps speed recovery from jet-lag. A session in the evening will delay your sleep cycle, a session in the morning will have an opposite effect.

IMPORTANT NOTES ON BRIGHT LIGHT THERAPY

The lamp comes with an electronic ballast which raises the light frequency up to 30 000 Hz, where flicker is not noticeable to the eye. After the first few sessions, mild headache or eye strain may occur. This is quite normal and such side effects usually disappear in a few days. If this does not happen, stop the therapy and consult your doctor before continuing. Not all people suffer from SAD due to lack of daylight. Exhaustion, sleep disorders and depression may have other causes. Bright light therapy is most effective when combined with regular physical exercise and a healthy diet.

Consult your doctor before starting bright light therapy if you suffer from

- eye illnesses (e.g. eye inflammation or hypersensitivity)
- mental disorders
- or if you use psychotropic drugs

When you start the bright light therapy, it is useful to take notes on the effects it has on you. This helps you determine which way the therapy works best for you. The need for bright light therapy as well as its duration and frequency vary individually. Use the therapy diary on the back of this booklet.

Congratulations for purchasing INNOLUX bright light therapy lamp. Before using the lamp, it is good to know what effects bright light has on your body and how to use bright light therapy effectively. Many people suffer from lethargy and depression during the dark winter season. This condition is known as SAD (Seasonal Affective Disorder). Scientific research has shown that as the amount of natural light decreases, the body may physiologically respond to these "unfavourable living conditions" by adjusting certain hormone levels and biological activities. SUB SAD or the Winter Blues is a milder form of SAD affecting a considerably larger number of people. The effects of these adjustments may include:

- Intense fatigue, exhaustion and various sleep disorders
- Increased appetite leading to weight gain
- General passivity and depression
- Decreased work efficiency and loss of vitality

On a sunny day the strength of light is approximately 50 000 – 100 000 lux, whereas indoors it is as little as 100-200 lux. An effective bright light therapy session requires the eyes to be exposed to at least 2 500 lux. Via the retina and the optic nerves the light affects the brain's light sensitive region (suprachiasmatic cellkernel), which adjusts the body's physiological cycles in accordance with the rhythm of day and night. The light reduces the amount of the sleep hormone melatonin, secreted by the pineal gland. Melatonin is a key factor in the regulation of alertness and sleep cycles, but also in the prevention of SAD. Bright light therapy helps you alleviate these problems and regain vitality. If the symptoms are severe, they should be treated under medical supervision.

Start bright light therapy before any symptoms of SAD appear, preferably soon after nights start drawing in.

Take a bright light session 3-5 times a week throughout the winter. If the symptoms are not relieved or they reappear, try taking daily sessions for 1-2 weeks. You may even begin with daily sessions or repeat this intensity whenever you feel you need it.

One session lasts for 0,5-2 hours. Longer than two hours does not intensify the effect but there's no harm of longer sessions. It is advisable to take the session between 6.00 a.m. and 6.00 p.m. but it is most effective in the morning between 6 and 8. A session in the evening may disturb sleep. If you tend to have evening sleepiness (common to elderly people), you might want to try taking the sessions earlier in the evening in order to stay up longer. You can also divide the session into several shorter periods during the same day if it suits your schedule better.

There is no need to look straight into the light source, but you should keep it within your visual field during the whole session. You can read, write or do other things in the meantime. The effectiveness of bright light therapy correlates with the distance between your face and the light source. Make sure that your face is not further from the lamp than the announced therapy distance on the model instructions. You can shorten the required session time by sitting closer to the lamp.

Bright light therapy also helps speed recovery from jet-lag. A session in the evening will delay your sleep cycle, a session in the morning will have an opposite effect.

IMPORTANT NOTES ON BRIGHT LIGHT THERAPY

The lamp comes with an electronic ballast which raises the light frequency up to 30 000 Hz, where flicker is not noticeable to the eye. After the first few sessions, mild headache or eye strain may occur. This is quite normal and such side effects usually disappear in a few days. If this does not happen, stop the therapy and consult your doctor before continuing. Not all people suffer from SAD due to lack of daylight. Exhaustion, sleep disorders and depression may have other causes. Bright light therapy is most effective when combined with regular physical

exercise and a healthy diet.

Consult your doctor before starting bright light therapy if you suffer from

- eye illnesses (e.g. eye inflammation or hypersensitivity)
- mental disorders
- or if you use psychotropic drugs

When you start the bright light therapy, it is useful to take notes on the effects it has on you. This helps you determine which way the therapy works best for you. The need for bright light therapy as well as its duration and frequency vary individually. Use the therapy diary on the back of this booklet.

1. Install the fluorescent tubes according to the model instruction. In some of the models the tubes are already installed.
2. The lamp is not designed for use in a bathroom, sauna or any other damp place, nor outdoors.
3. Connect the plug to a wall socket. The charge must be 220 - 240 V.
4. Use the switch in the electric cord to turn on the lamp. The light levels of the new fluorescent tubes may fluctuate during the first few hours, but this usually ceases after approximately ten hours of use.
5. Switch off the power after use.

IMPORTANT!

To get your new dimmable Bright light therapy lamp work at its best, have the lamps on on full power for the first 10 to 20 hours.

REPLACING THE FLUORESCENT TUBES

1. When replacing the tubes, check that the new ones have similar specifications to the originals.
2. Place the lamp on its back on a dry, clean and flat surface. Make sure that the plug is disconnected.
3. Remove the protective shield and change the fluorescent tubes according to the model instructions sheet.
4. Once the fluorescent tubes are firmly in place, check that they work before replacing the shade.
 - Please do not discard tubes with normal household refuse.
5. Dispose of the tubes appropriately so as not to endanger people or the environment.

MAINTENANCE OF THE LAMP

Keep the lamp clean and dust free. Always disconnect and allow tubes to cool down before cleaning. Use a damp cloth to clean its outer surfaces as well as the fluorescent tubes and inner surfaces. Do not use any solvents or detergents which may harm the surfaces. Do not let any water leak into the lamp.

IMPORTANT

- Do not let children play with the lamp
- If the electric cord is damaged, it must be replaced by an authorised electrician

TROUBLESHOOTING**1. THE BRIGHT LIGHT DEVICE IS NOT TURNING ON**

The fluorescent tubes may loosen from the lampholders during transport. In all of the models, tubes are not even pre-attached on their place. Check whether the fluorescent tubes are firmly attached on the lampholders.

2. THE LIGHT FLICKERS

If the light flickers, then check that the dimmer switch is at maximum level. Dimmable tubes should be burned on full power for the first 10-20 hours. After that the light should not anymore flicker when dimmed.

3. LIGHT SEEMS DIM?

If your bright light seems dim at first, it is due to the fact that fluorescent tubes should be burned about 100 hours before they reach full efficiency. Keep the light on for approximately 100 hours.

1. Install the fluorescent tubes according to the model instruction. In some of the models the tubes are already installed. 2. The lamp is not designed for use in a bathroom, sauna or any other damp place, nor outdoors. 3. Connect the plug to a wall socket. The charge must be 220 - 240 V. 4. Use the switch in the electric cord to turn on the lamp. The light levels of the new fluorescent tubes may fluctuate during the first few hours, but this usually ceases after approximately ten hours of use. 5. Switch off the power after use.

IMPORTANT! To get your new dimmable Bright light therapy lamp work at its best, have the lamps on on full power for the first 10 to 20 hours.

1. When replacing the tubes, check that the new ones have similar specifications to the originals. 2. Place the lamp on its back on a dry, clean and flat surface. Make sure that the plug is disconnected. 3. Remove the protective shield and change the fluorescent tubes according to the model instructions sheet. 4. Once the fluorescent tubes are firmly in place, check that they work before replacing the shade.

- Please do not discard tubes with normal household refuse. 5. Dispose of the tubes appropriately so as not to endanger people or the environment.

Keep the lamp clean and dust free. Always disconnect and allow tubes to cool down before cleaning. Use a damp cloth to clean its outer surfaces as well as the fluorescent tubes and inner surfaces. Do not use any solvents or detergents which may harm the surfaces. Do not let any water leak into the lamp.

IMPORTANT

- Do not let children play with the lamp
- If the electric cord is damaged, it must be replaced by an authorised electrician

1. THE BRIGHT LIGHT DEVICE IS NOT TURNING ON

The fluorescent tubes may loosen from the lampholders during transport. In all of the models, tubes are not even pre-attached on their place. Check whether the fluorescent tubes are firmly attached on the lampholders.

2. THE LIGHT FLICKERS

If the light flickers, then check that the dimmer switch is at maximum level. Dimmable tubes should be burned on full power for the first 10-20 hours. After that the light should not anymore flicker when dimmed.

3. LIGHT SEEMS DIM?

If your bright light seems dim at first, it is due to the fact that fluorescent tubes should be burned about 100 hours before they reach full efficiency. Keep the light on for approximately 100 hours.

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